



For the use of registered medical practitioners or hospitals or a laboratory



URINE CHEMISTRY ANALYZER

UroColor Series

UroMeter 720

UroMeter 120

UroMeter 120

Urine Chemistry Analyzer

GENERAL INFORMATION

- Auto calibration after startup
- Simple connectivity via USB communication port
- Ergonomic design enhances user convenience even for left-handed users
- Enhance work efficiency by easily switching between quick mode and normal mode
- Easy to input multiple IDs by keyboard, PC and barcode reader
- Easily switch to your language by user-defined method
- 4 ~ 11 parameter testing available
- Maximum 300 tests/hour, average 120 tests/hour
- Easy to learn and operate



VARIOUS INTERFACE PORTS

- RS232C, USB keyboard port, USB port
- Upload/download test results via bi-directional interface
- Handheld barcode reader for identifying specimen

WIDE LCD DISPLAY WINDOW

- Bright and clear color LCD
- 320*240 resolution
- Automatic switching to testing mode from standby mode

FAST RESULT PRINTING

- High-speed thermal printer
- Complete data about patient and test result

COMPACT SIZE, BUT HIGH CAPACITY

- Lightweight: 1.2 kg (2.65 lb)
- Compact size: 252 × 200 × 114 mm (W × D × H)

ERGONOMIC DESIGN

- Designed to streamline workflow
- Easy to use even for left-handed users

SPECIFICATIONS

- Power Adapter: DC 12V 3.33A (100–240V, 50/60 Hz)
- Size: 252 × 200 × 114 mm (W × D × H)
- Weight: 1.2 kg (2.65 lb)
- Memory: 100,000 results
- Print: Internal thermal printer
- LCD: 320*240 QVGA color
- RS232C Interface:
 - COM1: Communication with PC
 - COM2: Communication with barcode system

UroColor

Urine Chemistry Strip

TEST FOR BLOOD, BILIRUBIN, UROBILINOGEN, KETONE, PROTEIN, NITRITE, GLUCOSE, pH, SPECIFIC GRAVITY, LEUCOCYTES AND ASCORBIC ACID IN URINE.

GENERAL INFORMATION

- Clear positive and negative results
- Checks for interference in glucose, blood, bilirubin and nitrite by Vitamin C ingested
- Accurate results in about 60 seconds
- Shelf life: 24 months from the date of manufacturer
- Different color packages by parameters
- Fast result visually or instrumentally
- Packing size: 100's



UROCOLOR PARAMETERS

Cat No.	Product	Blood	Bilirubin	Urobilinogen	Ketone	Protein	Nitrite	Glucose	pH	Specific Gravity	Leucocyte	Ascorbic Acid
10UK02	UroColor 2					•		•				
10UK10	UroColor 10	•	•	•	•	•	•	•	•	•	•	

ORDERING INFORMATION

Product	Cat No.	Type	Pack Size
UroMeter 720	UM0720	Instrument	1 Unit
UroMeter 120	UM0120	Instrument	1 Unit
UroColor Series	Refer to table of UroColor parameters	Strip	100 Tests/Kit

FOR MORE INFORMATION CONTACT US:

Abbott Diagnostics Medical Private Limited.

No: 404, 4th Floor, BPTP Park Centra, NH-8, Opp. 32nd Milestone Gurugram-122 001, Haryana, India

Office No.: +91-124-4569000, E-mail: india.custsupport@abbott.com

© 2022 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. COL-10956-01 02/22

